

Feasibility and Efficacy of Three Methods of Zebra Mussel Larvae Detection in Lake Winnipeg


Sharon Lavigne, Mattias Johansson, Daniel Heath, Hugh MacIsaac

GLIER, University of Windsor



Dreissena spread

2013



Lake Winnipeg

- Manitoba
- Two basins
 - 2013 detected in south basin
 - 1 veliger reported in the north basin in October 2015
- Zebra mussels
 - 2013 detected in south basin
 - 1 veliger reported in the north basin in October 2015
- Low abundance throughout



Janusz (2015), personal communication

Project goals

- Detection of zebra mussels in Lake Winnipeg
- Comparison of methods of detection




Null Hypotheses

- Increased sampling intensity does not change probability of detection of NIS
- No difference in detection success in the south and north basin
- Different methods of detection would result in the same success rates

Methods

- July 2015
- Sampling sites
 - Grand Rapids
 - Dauphin River
 - Calder's Dock
 - Hnausa
 - Winnipeg Beach




Methods

- 100 plankton tows at each site
 - 64 μ m plankton nets
- Plankton collected and stored in 95% ethanol
- Filtered out plankton larger than 300 μ m
- Concentrated into 50mL




Analysis



Conventional Microscopy


- Cross-polarized light microscopy
 - Entire sample analyzed and enumerated
 - Distinct shape



Frischer *et al.* (2012); Johnson (1995)


Flow Cytometry and Image Analysis

- Flow Cytometer And Microscopy (FlowCAM)
 - Fitted with XPL (cross-polarized lenses)
 - Subsample (3/50 of filtered)
 - Images of each particle is captured and stored for analysis





Environmental DNA (eDNA)


- DNA detection
 - Subsample (1/50) of filtered plankton
 - DNA extraction
 - PCR amplified with species-specific primers




Microscopy: Abundance



Flow Rate




Increased Mortality




Horvath & Crane (2010)


Veliger Prevalence in Samples




Prevalence
(%)




Veliger Prevalence in Samples




Prevalence
(%)




Veliger Prevalence in Samples



Prevalence
(%)



Veliger Detection



Logistic Regression

Cost and time

	Microscopy	FlowCAM	eDNA
Start up cost (CAD\$)	10,441.60	81,499.73	5,835.00

Cost and time

	Microscopy	FlowCAM	eDNA
Start up cost (CAD\$)	10,441.60	81,499.73	5,835.00
Cost per sample (CAD\$)	7.90	6.78	1.87

Cost and time

	Microscopy	FlowCAM	eDNA
Start up cost (CAD\$)	10,441.60	81,499.73	5,835.00
Cost per sample (CAD\$)	7.90	6.78	1.87
Time (Hours)	198	167	24

Limitations

- Possible artifact of subsampling
- FlowCAM – volume restricted by time taken to process
- eDNA – volume restricted by equipment

Significance

- Changes in sampling strategy
- Detection methods
 - Microscopy mid-range in cost and longest time taken, but highest sensitivity
 - eDNA lowest cost and fastest - relatively high sensitivity

Acknowledgements

- Dr. Doug Haffner
- Dr. Aaron Fisk
- Dr Charles Ramcharan
- Colin van Overdijk
- MacIsaac Lab
- Heath Lab
- Joelle Pecz
- Jazel Gicole
- Taylor Haludek
- Ian MacIsaac
- Sarah-Jayne Collins



Multiple Stressors and Cumulative Effects in the Great Lakes:

An NSERC CREATE Program to Develop Innovative Solutions through International Training Partnerships